

# Optimization Study - Differentiating pancreatic cancer from individuals with concerning symptoms, including type II diabetes

Payam Delfani<sup>2</sup>, Anders Carlsson<sup>2</sup>, Thomas King<sup>2</sup>, Randall Brand<sup>3</sup>, Alexander Ney<sup>1</sup>, Stephen P Pereira<sup>1</sup>, Linda Dexlin Mellby<sup>2</sup>


1. University College London, UK 2. Immunovia, Lund, Sweden and Marlborough, USA 3. Department of Medicine, University of Pittsburgh, USA

## Objectives

IMMray™ PanCan-d Optimization Study aimed to evaluate how IMMray™ biomarker signature could separate patients with PDAC (stage I-IV) from individuals with various concerning symptomatic conditions not caused by PDAC, which mirrors the clinical setting encountered by healthcare professionals.

## Patients and Methods


In total, 923 serum samples were analyzed with IMMray™ discovery set up and a CA 19-9 ELISA. Patient samples from 136 PDAC (stage I-IV), 570 symptomatic individuals and 217 healthy controls were tested in a blinded manner. All PDACs were histologically confirmed. Based on one year follow up data non of the symptomatic controls developed PDAC. To minimize confounding and pre-analytical variables, all patient samples were collected and processed using the same standard operating procedures, stored at -80°C and tested within a year after collection. Data analysis for each group was performed using Support Vector Machine (SVM) algorithms. Data was divided into a training and test set, and test performance given as ROC AUC values, was then evaluated for the test set.




PDAC				Controls		
Stage I	Stage II	Stage III	Stage IV	Healthy controls	Symptomatic controls (without diabetes)	Diabetes controls
No. 20	34	21	61	217	480	90

## Results


### A. Symptomatic vs. PDAC




### B. Healthy vs. PDAC



### C. Diabetes vs. PDAC



### D. Controls vs. PDAC Stages I & II



**Fig 1.** In total, 923 individuals were analyzed. Combining IMMray™ biomarker signature with CA 19-9, the results from the test set showed ROC AUC values of 0.973 and 0.982 differentiating PDAC (stage I-IV) vs. symptomatic individuals, and healthy controls, respectively. PDAC (stage I-IV) could also be discriminated vs. diabetes type II controls with a ROC AUC value of 0.966.

**Fig 2.** PDAC Stages I & II could be discriminated from controls (symptomatic + healthy + diabetes), using IMMray™ biomarker signature and CA 19-9 ELISA. The result from the test set showed a ROC AUC value of 0.984 differentiating PDAC Stage I & II.

## Conclusions

- IMMray™ PanCan-d **Optimization Study** showed for the first time that IMMray™ biomarker signature together with CA 19-9 has the capacity to differentiate PDAC (stage I-IV) from symptomatic, non-PDAC individuals, including type II diabetes. This study paves the way for the next, where the IMMray™ PanCan-d commercial biomarker signature is selected and the commercial test model is built.
- Importantly, early stage I & II of PDAC was discriminated from controls with an unprecedented accuracy of 0.984.
- These findings need to be validated but have significant clinical implications for individuals attending primary and secondary care units with non-specific but concerning symptoms where PDAC may be suspected.

**Acknowledgements:** Immunovia would like to acknowledge Växjö Central Hospital, Department of Transfusion Medicine, Sweden, for providing freshly collected samples for the optimization study.

The IMMray™ PanCan-d optimization study was finalized 2019.

## References

1. Mellby LD, Nyberg AP, Johansen JS, et al. Serum Biomarker Signature-Based Liquid Biopsy for Diagnosis of Early-Stage Pancreatic Cancer. *Journal of Clinical Oncology*. 2018 Aug;36(28):2887-2894
2. Carl A. K. Borrebaeck. Precision diagnostics: moving towards pro marker signatures of clinical utility in cancer. *Nature Reviews Cancer*. 2017 Mar;17(3):199-204.